

Vision® 1327 Cart and Utensil Washer/Disinfector

Introduction This online training, for equipment technicians, will provide you with the knowledge to operate, maintain, service and troubleshoot the Vision® 1327 Cart and Utensil Washer/Disinfector.

Course Goal This course provides the student with the knowledge of how to operate the Vision® 1327 Cart and Utensil Washer/Disinfector and perform troubleshooting, preventive maintenance and adjustments. The student will learn the theory of operation, normal and service modes, functional descriptions of the major components and which components must be checked according to the Preventive Maintenance Guide.

Course Outcomes At the end of this course, the student will be able to:

- list the eight utilities required to operate the Vision® 1327
- locate the power on/off toggle switch
- locate the two pneumatic door cylinders and describe their purpose and function
- locate the two door gaskets and describe their purpose and function
- locate the two emergency stop cables and describe their purpose and function
- identify the six primary safety components
- identify the five components that prevent cross-contamination
- identify the guidelines for safe operation
- perform recommended checks for daily unit inspection prior to operation
- describe how to use the touch screen to access the Cycle mode and operate cycles
- describe how to use the touch screen to access supervisor mode to make changes
- locate the three RTDs and describe their purpose and function
- locate the three types of steam traps and describe their purpose and function
- locate the heat exchanger and describe its purpose and function
- locate the sump filter and describe its purpose and function
- locate the vacuum breaker and describe its purpose and function
- locate the safety valve and describe its purpose and function
- locate the water valves and describe their purpose and function
- locate the flow reducers and describe their purpose and function
- use the Preventive Maintenance Guide to identify and perform preventive maintenance procedures
- locate the Cycle Description Chart and describe its purpose
- list the schematic diagrams in the Reference Drawings section of the maintenance manuals
- identify two methods to use to contact STERIS technical support

Continued on next page

Vision® 1327 Cart and Utensil Washer/Disinfector, Continued

Why train online

You can decide the schedule and let your staff learn at their own pace. The training can be stopped as needed and picked up from where it is left off. Without ever leaving your facility, your staff can gain new skills, boost their confidence, and improve their productivity. An added benefit is that there are no travel costs.

Equipment Needed

Computer with internet access.

This online training is a self-contained course and does not include a maintenance manual as part of the cost of the course. If you would like to purchase a maintenance manual for this product, please contact STERIS Customer Support to place your order.

Target Audience

Biomedical Equipment Technicians and/or Hospital Maintenance Technicians

Duration

1 - 2 Hours

Assessment

End of course assessment required for completion.

Evaluations

End of course evaluation

Additional Information

<http://servicetraining.steris.com/>
