

STERIS® V-PRO® 60 Low Temperature Sterilization System

Introduction This online training is for equipment technicians who have previously trained on the STERIS® V-PRO® Series Low Temperature Sterilization Systems and are authorized, by their facility, to service the V-PRO® 60 Low Temperature Sterilization System. This online training will provide you with the differences between the V-PRO® 60 Low Temperature Sterilization System and the STERIS® V-PRO® 1, V-PRO® 1 Plus, and V-PRO® maX Series Low Temperature Sterilization Systems.

Course Goal This course provides an overview of how to operate the V-PRO® 60 Low Temperature Sterilization System. The student will learn the theory of operation, normal, supervisor, and service modes, functional descriptions of the major components and which components must be checked according to the Preventive Maintenance Schedule.

Course Outcomes At the end of this course, the student will be able to:

- identify the safety precautions to be taken when operating and servicing the V-PRO® 60 Low Temperature Sterilization System
- identify the four components of the front panel assembly
- describe how to operate the system in normal, supervisor, and service modes
- list the steps necessary to power on the control
- list the steps required to initiate a cycle
- list the steps required to enter supervisor mode
- list the steps required to enter service mode
- list the diagnostic and processing cycles
- identify the main components of the control and power box assemblies
- identify the pathway of VAPROX through the system components using the piping schematic
- list the six additional maintenance procedures to be performed during the annual inspection
- identify and describe the purpose and function of the draw tube assembly, injection manifold, injection cylinder and chamber solenoid valves
- access the STERIS® service websites for technical information and assistance

Why train online You can decide the schedule and let your staff learn at their own pace. The training can be stopped as needed and picked up from where it is left off. Without ever leaving your facility, your staff can gain new skills, boost their confidence, and improve their productivity. An added benefit is that there are no travel costs.

STERIS® V-PRO® 60 Low Temperature Sterilization System,

Continued

Equipment Needed

Computer with internet access.

This online training is a self-contained course and does not include a maintenance manual as part of the cost of the course. If you would like to purchase a maintenance manual for this product, please contact STERIS Customer Support to place your order.

Prerequisites

STERIS® V-PRO® 1, V-PRO® 1 Plus, or V-PRO® maX Low Temperature Sterilization Systems instructor led classroom training

Target Audience

Biomedical Equipment Technicians and/or Hospital Maintenance Technicians

Duration

1 – 2 Hours

Assessment

End of course assessment required for completion.

Evaluations

End of course evaluation

Additional Information

<http://servicetraining.steris.com/>
